

***There is only one
Thermo Scientific ROSS®***

Accept no Substitutes

**The pH electrode
for every application—
Thermo Scientific
ROSS Electrodes**

Experience the accuracy
and trouble-free results
from using the ideal electrode
suited for your application.



***Purchase a
Thermo Scientific
ROSS Electrode
and receive a FREE
Buffer Calibration kit!***

FREE pH Buffer Calibration Kit (810199):

475 mL bottles of 4, 7 and 10 Buffers, Cleaning Solution and ROSS Storage Solution

Electrode Application Guide

APPLICATION	Cat No	Description	Advantage
Cheese and Meats Agar plates	8135BN	ROSS Flat Surface, epoxy-body electrode	Flat surface for surface contact, and no silver in reference to interfere with possible sulfates or proteins.
Creams, Gels and Emulsions	8172BNWP	ROSS Sure-Flow glass-body electrode	Free-flowing junction assures the most stable, drift-free measurements. Easy to clean junction never clogs.
Drinking water General Purpose	8104BN	ROSS combination glass-body with rugged bulb	The features and accuracy of a ROSS with a durable sensing bulb.
Low ionic strength: Boiler water, Well water, Distilled water	8157BNUMD	ROSS Ultra pH/ATC epoxy-body electrode	The features and accuracy of a ROSS with the convenience of a triode.
Pharmaceutical	8102BNUWP	ROSS Ultra combination glass-body	Superior stability, accuracy and stability with best temperature performance. Tris, protein and sulfide compatible.
High viscosity: Slurries and Sludges	8165BNWP	ROSS Sure-Flow epoxy-body electrode	Free-flowing junction assures the most stable, drift-free measurements. Easy to clean junction never clogs.
Wastewater Seawater	8172BNWP	ROSS Sure-Flow glass-body electrode	Free-flowing junction assures the most stable, drift-free measurements. Easy to clean junction never clogs.

Note: We offer a range of electrodes per application—Call our technical support line for specific application challenges at 1-800-225-1480 or check our application notes at www.thermo.com/water.

It's easy to redeem:

Simply place your order, complete the redemption form and fax along with your invoices or packing slips to 978.232.1042.

Your Name:		
Company:		
Street Address:		
City:	State:	Zip Code:
E-mail Address:		
Telephone:		

Limited time only, while supplies last, offer commencing May 1, 2010 and is valid in North America. Product will be drop shipped from the factory to the end user. This offer may not be combined with any other coupon, discount or package. All prices are in U.S. dollars. This offer is void where prohibited by law or company policy. Offer expires December 31, 2010. Please allow 4–6 weeks for delivery of promotional product.



P-ROSS-E 0310 RevA

Canadawide
SCIENTIFIC

Thermo
SCIENTIFIC